



**CUSTOM MADE POST-LIPOSUCTION COMPRESSION GARMENTS**

**MEASUREMENT CHART**

**MEASUREMENTS  
(IN CENTIMETERS)  
CIRCUMFERENCE**

**INSTRUCTIONS**  
(PLEASE READ THE FOLLOWING INSTRUCTIONS CAREFULLY BEFORE FILLING IN THE MEASUREMENTS.)

1. Fill in all the measurements in CENTIMETRES.
2. Select your design and fill in the required measurements only. Shaded block means that the measurement is not required.
3. Take the measurements on Bare Skin ( i.e. without the clothes.)
4. For length measurements keep the limbs straight.
5. Use ordinary measuring tape.

- For the convenience of adjustments Post Liposuction Compression Garments are made with Velcro opening on both lateral sides.

Only subsequent garments can be made with zipper.

ANY ADDITIONAL INFORMATION

---




---



---



---



REQUIRED MEASUREMENTS	(116)	(117)	(118)	(119)
(A) TOP OF GARMENT				
(B) WAIST ON NAVEL				
(C) MAXIMUM HIP				
(D) GLUTEAL FOLD				
(E) MIDDLE OF THIGH				
(F) JUST ABOVE KNEE				
(G) AT KNEE				
(H) JUST BELOW KNEE				
(I) MAXIMUM CALF				
(J) END OF GARMENT				

**LENGTH**

(1) A To B				
(2) B To D				
(3) D TO F				
(4) D TO G				
(5) G TO J				

DESIGN CHOICES	YES
OPEN PUBIS	
PANTY LT. OPEN	
PANTY RT. OPEN	
SHOULDER STRAPS	

MAIL THIS FORM & D.D. OF THE APPROPRIATE AMOUNT TO :

**NORMA D.N.D. PRODUCTS PVT. LTD.**

NORMA Complex, D.D.A. Market, J-Block, Vikas Puri,  
N.Delhi - 110 018, PH: **2854-1111**(5 Lines),

**98182 99999**(5 Lines), FAX : 011-28541122

E-mail: normadnd@vsnl.com Website: www.normadnd.com